

# SYNCopen™ for CT motion™

Interface between the contrast media injector and CT scanner

## What are the advantages of SYNCopen?

- **Real-time data exchange** via the interface between the contrast media injector and CT scanner
- Quality assurance through the **automatic sending of injection parameters** from the injector to the CT scanner
- **Improved workflow** thanks to a single-button start



## Get more out of your CT motion

Function	CT motion without SYNCopen	CT motion with SYNCopen
Synchronous start & stop of injector and scanner	-	■
Provision and display of injection protocols on the scanner	-	■
Transmission of injection parameters to the RIS/PACS system	-	■
Optimized image quality thanks to synchronous start	-	■
One-time entry of injection protocols saves time and administrative effort	-	■

SYNCopen supports  
CANopen Class 4