

# Lightweight and high-performance

Heating cabinet for contrast media heating

## Compact heating cabinet for flexible use

- Small and lightweight for CT use
- Heating of viscous contrast media for higher flow rates
- Ready for immediate use thanks to easy handling
- Optimized patient comfort through the use of temperature-controlled contrast media

## Device Features

- Electronic temperature control
- Alarm in the event of a malfunction
- Rustproof stainless steel housing
- Glass door with viewing window, lockable
- Temperature range up from 30 °C to 70 °C
- 2 stainless steel shelves, variable applications



## Technical Data

<b>Dimensions (W × D × H)</b>	560 mm × 370 mm × 380 mm
<b>Weight</b>	10.2 kg incl. perforated plates
<b>Capacity</b>	20 l
<b>Power Connection</b>	230 V
<b>Frequency</b>	50/60 Hz
<b>Power input</b>	100 W