

# uCentum™

comprehensive posterior system

# uBase™

iliosacral fixation system



Complete  
system for open  
and percutaneous  
applications


**Ulrich**  
medical

*dedicated to you.*

# uCentum™

## System


- Posterior stabilization of the thoracic to lumbosacral spine
- Complete system for open and percutaneous applications
- Ø 5.5 mm and Ø. 6.0 mm rods
- Polyaxial and monoaxial screws
- Pre-fixation screws and long-arm screws
- Can be combined with all the screws from the uBase system
- Use in combination with the Brainlab navigation software possible (both open and percutaneous)
- Extensive navigation module for preparation of the surgical procedure and screw insertion
- Combination with the cervicothoracic screw-rod system neon<sup>3</sup>
- Implants available sterile



Cannulated and perforated screws (except Ø 4.5 mm)




Sleeve for percutaneous application



Offset connector

Low-profile screw head

High flexibility due to various  
connector lengths and additional  
connection options



# uBase™

## System

- Stabilization of the sacrum and ilium in connection with uCentum
- Purely sacral, sacroiliac or purely iliac placement
- Polyaxial screw head, pre-angled medial/lateral or cranial/caudal
- Additionally available: fixed-angle screw
- All screws cannulated
- Implants available sterile

## Advantages


- Improved safety thanks to cannulated screws and extensive navigation module for preparation of the surgical procedure and screw insertion
- 48° polyaxiality
- Percutaneous application (polyaxial screws)
- Screw head with low profile




Iliosacral screw, polyaxial



Purely sacral screw placement



Sacroiliac screw placement




Purely iliac screw placement



# uCentum™

## Advantages

- Improved safety thanks to cannulated screws and extensive navigation module for preparation of the surgical procedure and screw insertion
- Cement augmentation through the screw thanks to lateral perforations
- Rods with and crosslinks for Ø 5.5 and Ø 6.0 mm
- Combination with the cervicothoracic screw-rod system neon<sup>3</sup>
- New reduction instrumentation



# Components

## Implants

Screw, polyaxial, Ø 4.5 mm, length 25 mm – 50 mm  
 Screw, polyaxial, Ø 5.5 mm, length 25 mm – 55 mm, perforated  
 Screw, polyaxial, Ø 6.5 mm, length 30 mm – 60 mm, perforated  
 Screw, polyaxial, Ø 7.5 mm, length 35 mm – 60 mm, perforated  
 Screw, polyaxial, Ø 8.5 mm, length 35 mm – 60 mm, perforated  
 Screw, polyaxial, Ø 10 mm, length 45 mm – 60 mm, perforated  
 Screw, polyaxial, sterile, Ø 11.5 mm, length 45 mm – 55 mm, perforated\*

## Product no.

CS 3802-045-025 – CS 3802-045-050  
 CS 3802-055-025 – CS 3802-055-055  
 CS 3802-065-030 – CS 3802-065-060  
 CS 3802-075-035 – CS 3802-075-060  
 CS 3802-085-035 – CS 3802-085-060  
 CS 3802-100-045 – CS 3802-100-060  
 CS 3802-115-045-S – CS 3802-115-055-S

Screw, monoaxial, Ø 4.5 mm, length 25 mm – 50 mm  
 Screw, monoaxial, Ø 5.5 mm, length 25 mm – 55 mm, perforated  
 Screw, monoaxial, Ø 6.5 mm, length 30 mm – 60 mm, perforated  
 Screw, monoaxial, Ø 7.5 mm, length 35 mm – 60 mm, perforated  
 Screw, monoaxial, Ø 8.5 mm, length 35 mm – 60 mm, perforated  
 Screw, monoaxial, Ø 10 mm, length 45 mm – 60 mm, perforated

CS 3803-045-025 – CS 3803-045-050  
 CS 3803-055-025 – CS 3803-055-055  
 CS 3803-065-030 – CS 3803-065-060  
 CS 3803-075-035 – CS 3803-075-060  
 CS 3803-085-035 – CS 3803-085-060  
 CS 3803-100-045 – CS 3803-100-060

Screw, pre-fixation, Ø 4.5 mm, length 25 mm – 50 mm  
 Screw, pre-fixation, Ø 5.5 mm, length 25 mm – 55 mm, perforated  
 Screw, pre-fixation, Ø 6.5 mm, length 30 mm – 60 mm, perforated  
 Screw, pre-fixation, Ø 7.5 mm, length 35 mm – 60 mm, perforated  
 Screw, pre-fixation, Ø 8.5 mm, length 35 mm – 60 mm, perforated  
 Screw, pre-fixation, Ø 10 mm, length 45 mm – 60 mm, perforated

CS 3804-045-025 – CS 3804-045-050  
 CS 3804-055-025 – CS 3804-055-055  
 CS 3804-065-030 – CS 3804-065-060  
 CS 3804-075-035 – CS 3804-075-060  
 CS 3804-085-035 – CS 3804-085-060  
 CS 3804-100-045 – CS 3804-100-060

Long-arm screw, pre-fixation, Ø 5.5 mm, length 25 mm – 55 mm, perforated  
 Long-arm screw, pre-fixation, Ø 6.5 mm, length 30 mm – 60 mm, perforated  
 Long-arm screw, pre-fixation, Ø 7.5 mm, length 35 mm – 60 mm, perforated  
 Long-arm screw, pre-fixation, Ø 8.5 mm, length 35 mm – 60 mm, perforated

CS 3805-055-025 – CS 3805-055-055  
 CS 3805-065-030 – CS 3805-065-060  
 CS 3805-075-035 – CS 3805-075-060  
 CS 3805-085-035 – CS 3805-085-060

## Combination with uBase

Iliosacral screw, Ø 7.5 mm, length 40 mm – 80 mm  
 Iliosacral screw, Ø 8.5 mm, length 40 mm – 120 mm  
 Iliosacral screw, Ø 10.0 mm, length 80 mm – 140 mm

CS 7170-075-040 – CS 7170-075-080  
 CS 7170-085-040 – CS 7170-085-120  
 CS 7170-100-080 – CS 7170-100-140

Iliosacral screw, polyaxial, Ø 7.5 mm, length 40 mm – 80 mm, sterile\*  
 Iliosacral screw, polyaxial, Ø 8.5 mm, length 40 mm – 120 mm, sterile\*  
 Iliosacral screw, polyaxial, Ø 10.0 mm, length 80 mm – 140 mm, sterile\*

CS 7172-075-040-S – CS 7172-075-080-S  
 CS 7172-085-040-S – CS 7172-085-120-S  
 CS 7172-100-080-S – CS 7172-100-140-S

Iliac screw, polyaxial, sterile, Ø 8.5 mm, length 80 mm – 120 mm\*  
 Iliac screw, polyaxial, sterile, Ø 10.0 mm, length 80 mm – 120 mm\*

CS 7173-085-080-S – CS 7173-085-120-S  
 CS 7173-100-080-S – CS 7173-100-120-S

\*These products are only available sterile

## Information on sterile packed implants

To order sterile packed implants, add "S" to the end of the product number (example: CS 3802-045-025-S).



**ulrich GmbH & Co. KG**  
 Buchbrunnenweg 12  
 89081 Ulm, Germany  
 T +49 731 9654-0  
 spine@ulrichmedical.com  
 www.ulrichmedical.com

**ulrich medical France SAS**  
 10, Allée Aristide Maillol  
 31770 Colomiers, France  
 T +33 5 34 50 91 02  
 info@ulrichmedical.fr  
 www.ulrichmedical.fr

