

# ADD™

anterior distraction device




**Ulrich**  
medical

*dedicated to you.*

# System

- Distractable vertebral body replacement for the anterior cervical and cervico-thoracic spine
- Smoothly distractable for individual anatomical adaptation
- Large filling volume for bone/bone substitute material
- Clear instrument set and easy handling
- All implants also available in individual sterile packaging





## Advantages

- High product variability for individual, anatomical adaption
- Smooth, in situ distraction and compression
- Large filling volume



# Components

Implants	Product number
ADD, Ø 12 mm, angle 0°, height 10 – 13 mm	CS 2250-12-10
ADD, Ø 12 mm, angle 0°, height 12 – 17 mm	CS 2250-12-12
ADD, Ø 12 mm, angle 0°, height 16 – 25 mm	CS 2250-12-16
ADD, Ø 12 mm, angle 6°, height 16 – 25 mm	CS 2250-12-166
ADD, Ø 12 mm, angle 0°, height 24 – 40 mm	CS 2250-12-24
ADD, Ø 12 mm, angle 6°, height 24 – 40 mm	CS 2250-12-246
ADD, Ø 12 mm, angle 0°, height 39 – 65 mm	CS 2250-12-39
ADD, Ø 12 mm, angle 6°, height 39 – 65 mm	CS 2250-12-396
ADD, Ø 14 mm, angle 0°, height 10 – 13 mm	CS 2250-14-10
ADD, Ø 14 mm, angle 0°, height 12 – 17 mm	CS 2250-14-12
ADD, Ø 14 mm, angle 0°, height 16 – 25 mm	CS 2250-14-16
ADD, Ø 14 mm, angle 6°, height 16 – 25 mm	CS 2250-14-166
ADD, Ø 14 mm, angle 0°, height 24 – 40 mm	CS 2250-14-24
ADD, Ø 14 mm, angle 6°, height 24 – 40 mm	CS 2250-14-246
ADD, Ø 14 mm, angle 0°, height 39 – 65 mm	CS 2250-14-39
ADD, Ø 14 mm, angle 6°, height 39 – 65 mm	CS 2250-14-396
ADD, Ø 16 mm, angle 0°, height 10 – 13 mm	CS 2250-16-10
ADD, Ø 16 mm, angle 0°, height 12 – 17 mm	CS 2250-16-12
ADD, Ø 16 mm, angle 0°, height 16 – 25 mm	CS 2250-16-16
ADD, Ø 16 mm, angle 6°, height 16 – 25 mm	CS 2250-16-166
ADD, Ø 16 mm, angle 0°, height 24 – 40 mm	CS 2250-16-24
ADD, Ø 16 mm, angle 6°, height 24 – 40 mm	CS 2250-16-246
ADD, Ø 16 mm, angle 0°, height 39 – 65 mm	CS 2250-16-39
Locking screw	CS 2259

## Information on sterile packed implants

To order sterile packed implants, add "S" to the end of the product number (example: CS 2250-12-10-S).



**ulrich GmbH & Co. KG**  
 Buchbrunnenweg 12  
 89081 Ulm, Germany  
 T +49 731 9654-0  
 spine@ulrichmedical.com  
 www.ulrichmedical.com

**ulrich medical France SAS**  
 10, Allée Aristide Maillol  
 31770 Colomiers, France  
 T +33 5 34 50 91 02  
 info@ulrichmedical.fr  
 www.ulrichmedical.fr