

ADDplus™

anterior distraction device plus





Ulrich
medical

dedicated to you.

System


- Distractable vertebral body replacement with solid plate for the anterior cervical and cervicothoracic spine
- Increasing angles with increasing implant size (0° - 18°)
- Optionally cancellous bone screws or expandable screws
- Large filling volume for bone/ bone substitute material
- All implants also available in individual sterile packaging





Advantages

- Vertebral body substitute and plate, in a single system
- Optimized implant anchoring by expandable screw
- Smoothly distractable for individual anatomical adaptation
- Large filling volume



Components

Implants	Product number
ADDplus, Ø 12 mm, angle 0°, height 13 – 18 mm	CS 2253-12-13
ADDplus, Ø 12 mm, angle 6°, height 17 – 26 mm	CS 2253-12-17
ADDplus, Ø 12 mm, angle 12°, height 25 – 41 mm	CS 2253-12-25
ADDplus, Ø 12 mm, angle 18°, height 40 – 65 mm	CS 2253-12-40
osmium screw incl. bolt, gold, Ø 5.0 mm, length 14 mm	CS 1300-14T
osmium screw incl. bolt, green, Ø 5.0 mm, length 16 mm	CS 1300-16T
osmium screw incl. bolt, blue, Ø 5.0 mm, length 18 mm	CS 1300-18T
osmium screw incl. bolt, grey, Ø 4.0 mm, length 12 mm	CS 1303-12
osmium screw incl. bolt, gold, Ø 4.0 mm, length 14 mm	CS 1303-14
osmium screw incl. bolt, green, Ø 4.0 mm, length 16 mm	CS 1303-16
osmium screw incl. bolt, blue, Ø 4.0 mm, length 18 mm	CS 1303-18
Cancellous bone screw, Ø 4.0 mm, length 12 mm	CS 1301-12T
Cancellous bone screw, Ø 4.0 mm, length 14 mm	CS 1301-14T
Cancellous bone screw, Ø 4.0 mm, length 16 mm	CS 1301-16T
Cancellous bone screw, Ø 4.0 mm, length 18 mm	CS 1301-18T
Cancellous bone screw, Ø 4.0 mm, length 20 mm	CS 1301-20T
Cancellous bone screw, Ø 4.0 mm, length 22 mm	CS 1301-22T
Cancellous bone screw, Ø 4.0 mm, length 24 mm	CS 1301-24T
Cancellous bone screw, Ø 4.0 mm, length 26 mm	CS 1301-26T
Locking screw	CS 2259

Information on sterile packed implants

To order sterile packed implants, add "S" to the end of the product number (example: CS 2253-12-13-S).



ulrich GmbH & Co. KG
Buchbrunnenweg 12
89081 Ulm, Germany
T +49 731 9654-0
spine@ulrichmedical.com
www.ulrichmedical.com

ulrich medical France SAS
10, Allée Aristide Maillol
31770 Colomiers, France
T +33 5 34 50 91 02
info@ulrichmedical.fr
www.ulrichmedical.fr